



PLPP3 Antibody, HRP conjugated

Product Code	CSB-PA018415LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O14495
Immunogen	Recombinant Human Phospholipid phosphatase 3 protein (144-193AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Phospholipid phosphatase 3 (EC 3.1.3.4) (Lipid phosphate phosphohydrolase 3) (PAP2-beta) (Phosphatidate phosphohydrolase type 2b) (Phosphatidic acid phosphatase 2b) (PAP-2b) (PAP2b) (Vascular endothelial growth factor and type I collagen-inducible protein) (VCIP), PLPP3, LPP3 PPAP2B
Immunogen Species	Homo sapiens (Human)
Research Area	Cardiovascular
Target Names	PLPP3
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.