



PPIAL4G Antibody, HRP conjugated

| | |
|----------------------------|--|
| Product Code | CSB-PA018453LB01HU |
| Storage | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze. |
| Uniprot No. | P0DN37 |
| Immunogen | Recombinant Human Peptidyl-prolyl cis-trans isomerase A-like 4G protein (1-164AA) |
| Raised In | Rabbit |
| Species Reactivity | Human |
| Tested Applications | ELISA |
| Form | Liquid |
| Conjugate | HRP |
| Storage Buffer | Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4 |
| Purification Method | >95%, Protein G purified |
| Isotype | IgG |
| Clonality | Polyclonal |
| Alias | Peptidyl-prolyl cis-trans isomerase A-like 4G (PPIase A-like 4G) (EC 5.2.1.8) (Peptidylprolyl cis-trans isomerase A-like 4), PPIAL4G |
| Immunogen Species | Homo sapiens (Human) |
| Research Area | Others |
| Target Names | PPIAL4G |
| Usage | For Research Use Only. Not for use in diagnostic or therapeutic procedures. |