



# PPIAL4G Antibody, FITC conjugated

<b>Product Code</b>	CSB-PA018453LC01HU
<b>Storage</b>	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
<b>Uniprot No.</b>	P0DN37
<b>Immunogen</b>	Recombinant Human Peptidyl-prolyl cis-trans isomerase A-like 4G protein (1-164AA)
<b>Raised In</b>	Rabbit
<b>Species Reactivity</b>	Human
<b>Form</b>	Liquid
<b>Conjugate</b>	FITC
<b>Storage Buffer</b>	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
<b>Purification Method</b>	>95%, Protein G purified
<b>Isotype</b>	IgG
<b>Clonality</b>	Polyclonal
<b>Alias</b>	Peptidyl-prolyl cis-trans isomerase A-like 4G (PPIase A-like 4G) (EC 5.2.1.8) (Peptidylprolyl cis-trans isomerase A-like 4), PPIAL4G
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Research Area</b>	Others
<b>Target Names</b>	PPIAL4G
<b>Usage</b>	For Research Use Only. Not for use in diagnostic or therapeutic procedures.