



PPIAL4G Antibody, Biotin conjugated

Product Code	CSB-PA018453LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P0DN37
Immunogen	Recombinant Human Peptidyl-prolyl cis-trans isomerase A-like 4G protein (1-164AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Peptidyl-prolyl cis-trans isomerase A-like 4G (PPIase A-like 4G) (EC 5.2.1.8) (Peptidylprolyl cis-trans isomerase A-like 4), PPIAL4G
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	PPIAL4G
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.