




PPP2CA Antibody

Product Code	CSB-PA018559LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P67775
Immunogen	Recombinant Human Serine/threonine-protein phosphatase 2A catalytic subunit alpha isoform protein (1-309AA)
Raised In	Rabbit
Species Reactivity	Human, Mouse
Tested Applications	ELISA, WB, IHC, IF; Recommended dilution: WB:1:500-1:5000, IHC:1:20-1:200, IF:1:50-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Serine/threonine-protein phosphatase 2A catalytic subunit alpha isoform (PP2A-alpha) (EC 3.1.3.16) (Replication protein C) (RP-C), PPP2CA
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	PPP2CA

Image



Western Blot

Positive WB detected in: Mouse liver tissue,
Mouse kidney tissue


All lanes: PPP2CA antibody at 2.5µg/ml

Secondary


Goat polyclonal to rabbit IgG at 1/50000 dilution

Predicted band size: 36, 30 kDa

Observed band size: 36, 30 kDa



Immunohistochemistry of paraffin-embedded human liver cancer using CSB-PA018559LA01HU at dilution of 1:100



Immunofluorescent analysis of HeLa cells using CSB-PA018559LA01HU at dilution of 1:100 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.