




NPY4R Antibody

Product Code	CSB-PA018593LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P50391
Immunogen	Recombinant Human Neuropeptide Y receptor type 4 protein (328-375AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC; Recommended dilution: IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Neuropeptide Y receptor type 4 (NPY4-R) (Pancreatic polypeptide receptor 1) (PP1), NPY4R, PPYR1
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	NPY4R

Image



Immunohistochemistry of paraffin-embedded human breast cancer using CSB-PA018593LA01HU at dilution of 1:100

Usage For Research Use Only. Not for use in diagnostic or therapeutic procedures.