



NPY4R Antibody, HRP conjugated

Product Code	CSB-PA018593LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P50391
Immunogen	Recombinant Human Neuropeptide Y receptor type 4 protein (328-375AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Neuropeptide Y receptor type 4 (NPY4-R) (Pancreatic polypeptide receptor 1) (PP1), NPY4R, PPYR1
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	NPY4R
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.