



PRAMEF5 Antibody

Product Code	CSB-PA018622LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q5TYX0
Immunogen	Recombinant Human PRAME family member 5 protein (1-476AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	PRAME family member 5, PRAMEF5, PRAMEF23 PRAMEF5L
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	PRAMEF5
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.