




Prf1 Antibody

Product Code	CSB-PA018668LA01MO
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P10820
Immunogen	Recombinant Mouse Perforin-1 protein (21-554AA)
Raised In	Rabbit
Species Reactivity	Mouse
Tested Applications	ELISA, WB; Recommended dilution: WB:1:500-1:2000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	Antigen affinity purification
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Mus musculus (Mouse)
Research Area	Immunology
Target Names	Prf1

Image



Western Blot
 Positive WB detected in: CTLL2 whole cell lysate(30µg)
 All lanes: PRF1 antibody at 1:1000
 Secondary
 Goat polyclonal to rabbit IgG at 1/40000 dilution
 Predicted band size: 63 kDa
 Observed band size: 70 kDa
 Exposure time: 120s

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.