



Prf1 Antibody, HRP conjugated

Product Code	CSB-PA018668LB01MO
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P10820
Immunogen	Recombinant Mouse Perforin-1 protein (21-554AA)
Raised In	Rabbit
Species Reactivity	Mouse
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Perforin-1 (P1) (Cytolysin) (Lymphocyte pore-forming protein), Prf1, Pfp
Immunogen Species	Mus musculus (Mouse)
Research Area	Immunology
Target Names	Prf1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.