



Prnp Antibody, HRP conjugated

Product Code	CSB-PA018739LB01RA
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P13852
Immunogen	Recombinant Rat Major prion protein (29-231AA)
Raised In	Rabbit
Species Reactivity	Rat
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Major prion protein (PrP) (CD antigen CD230), Prnp, Prn Prp
Immunogen Species	Rattus norvegicus (Rat)
Research Area	Others
Target Names	Prnp