



Prnp Antibody, Biotin conjugated

Product Code	CSB-PA018739LD01RA
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P13852
Immunogen	Recombinant Rat Major prion protein (29-231AA)
Raised In	Rabbit
Species Reactivity	Rat
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Major prion protein (PrP) (CD antigen CD230), Prnp, Prn Prp
Immunogen Species	Rattus norvegicus (Rat)
Research Area	Others
Target Names	Prnp
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.