



PROKR1 Antibody, HRP conjugated

Product Code	CSB-PA018748LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8TCW9
Immunogen	Recombinant Human Prokineticin receptor 1 protein (1-62AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Prokineticin receptor 1 (PK-R1) (G-protein coupled receptor 73) (G-protein coupled receptor ZAQ) (GPR73a), PROKR1, GPR73 PKR1
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	PROKR1