



PROKR1 Antibody, FITC conjugated

Product Code	CSB-PA018748LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8TCW9
Immunogen	Recombinant Human Prokineticin receptor 1 protein (1-62AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Prokineticin receptor 1 (PK-R1) (G-protein coupled receptor 73) (G-protein coupled receptor ZAQ) (GPR73a), PROKR1, GPR73 PKR1
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	PROKR1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.