



Prss1 Antibody

Product Code	CSB-PA018811LA01RA
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P00762
Immunogen	Recombinant Rat Anionic trypsin-1 protein (24-246AA)
Raised In	Rabbit
Species Reactivity	Rat
Tested Applications	ELISA
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Anionic trypsin-1 (EC 3.4.21.4) (Anionic trypsin I) (Pretrypsinogen I) (Serine protease 1), Prss1, Try1
Immunogen Species	Rattus norvegicus (Rat)
Research Area	Cell Biology
Target Names	Prss1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.