



Prss1 Antibody, Biotin conjugated

Product Code	CSB-PA018811LD01RA
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P00762
Immunogen	Recombinant Rat Anionic trypsin-1 protein (24-246AA)
Raised In	Rabbit
Species Reactivity	Rat
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Anionic trypsin-1 (EC 3.4.21.4) (Anionic trypsin I) (Pretrypsinogen I) (Serine protease 1), Prss1, Try1
Immunogen Species	Rattus norvegicus (Rat)
Research Area	Cell Biology
Target Names	Prss1