



PRY Antibody, HRP conjugated

Product Code	CSB-PA018834LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O14603
Immunogen	Recombinant Human PTPN13-like protein, Y-linked protein (13-147AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	PTPN13-like protein, Y-linked (Testis-specific PTP-BL-related Y protein), PRY; PRY2; PRYP3; PRYP4, PRY1 PTPN13LY; PTPN13LY2
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	PRY
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.