



PSG6 Antibody, Biotin conjugated

Product Code	CSB-PA018856LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q00889
Immunogen	Recombinant Human Pregnancy-specific beta-1-glycoprotein 6 protein (35-300AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Pregnancy-specific beta-1-glycoprotein 6 (PS-beta-G-6) (PSBG-6) (Pregnancy-specific glycoprotein 6) (Pregnancy-specific beta-1-glycoprotein 10) (PS-beta-G-10) (PSBG-10) (Pregnancy-specific glycoprotein 10) (Pregnancy-specific beta-1-glycoprotein 12) (PS-beta-G-12) (PSBG-12) (Pregnancy-specific glycoprotein 12), PSG6, CGM3 PSG10 PSG12 PSGGB
Immunogen Species	Homo sapiens (Human)
Research Area	Tags & Cell Markers
Target Names	PSG6
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.