



PSMC3 Antibody, Biotin conjugated

Product Code	CSB-PA018891LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P17980
Immunogen	Recombinant Human 26S proteasome regulatory subunit 6A protein (1-300AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	26S proteasome regulatory subunit 6A (26S proteasome AAA-ATPase subunit RPT5) (Proteasome 26S subunit ATPase 3) (Proteasome subunit P50) (Tat-binding protein 1) (TBP-1), PSMC3, TBP1
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	PSMC3
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.