



QDPR Antibody, HRP conjugated

Product Code	CSB-PA019133LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P09417
Immunogen	Recombinant Human Dihydropteridine reductase protein (2-244AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Dihydropteridine reductase (EC 1.5.1.34) (HDHPR) (Quinoid dihydropteridine reductase) (Short chain dehydrogenase/reductase family 33C member 1), QDPR, DHPR SDR33C1
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	QDPR