



QDPR Antibody, HRP conjugated

Product Code	CSB-PA019133LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P09417
Immunogen	Recombinant Human Dihydropteridine reductase protein (2-244AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Dihydropteridine reductase (EC 1.5.1.34) (HDHPR) (Quinoid dihydropteridine reductase) (Short chain dehydrogenase/reductase family 33C member 1), QDPR, DHPR SDR33C1
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	QDPR
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.