



RANBP2 Antibody, Biotin conjugated

Product Code	CSB-PA019311LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P49792
Immunogen	Recombinant Human E3 SUMO-protein ligase RanBP2 protein (2601-2802AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	E3 SUMO-protein ligase RanBP2 (EC 2.3.2.-) (358 kDa nucleoporin) (E3 SUMO-protein transferase RanBP2) (Nuclear pore complex protein Nup358) (Nucleoporin Nup358) (Ran-binding protein 2) (RanBP2) (p270), RANBP2, NUP358
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	RANBP2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.