



RASA1 Antibody, FITC conjugated

Product Code	CSB-PA019346LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P20936
Immunogen	Recombinant Human Ras GTPase-activating protein 1 protein (201-500AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Ras GTPase-activating protein 1 (GAP) (GTPase-activating protein) (RasGAP) (Ras p21 protein activator) (p120GAP), RASA1, GAP RASA
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	RASA1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.