



# RASA1 Antibody, Biotin conjugated

Product Code	CSB-PA019346LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P20936
Immunogen	Recombinant Human Ras GTPase-activating protein 1 protein (201-500AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Ras GTPase-activating protein 1 (GAP) (GTPase-activating protein) (RasGAP) (Ras p21 protein activator) (p120GAP), RASA1, GAP RASA
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	RASA1