




RBL1 Antibody

Product Code	CSB-PA019398LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P28749
Immunogen	Recombinant Human Retinoblastoma-like protein 1 protein (201-500AA)
Raised In	Rabbit
Species Reactivity	Human, Mouse
Tested Applications	ELISA, WB; Recommended dilution: WB:1:2000-1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Retinoblastoma-like protein 1 (107 kDa retinoblastoma-associated protein) (p107) (pRb1), RBL1
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	RBL1

Image



Western Blot

Positive WB detected in: Mouse skeletal muscle tissue

All lanes: RBL1 antibody at 3µg/ml

Secondary

Goat polyclonal to rabbit IgG at 1/50000 dilution

Predicted band size: 121, 115 kDa

Observed band size: 121 kDa