



RBL1 Antibody, Biotin conjugated

Product Code	CSB-PA019398LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P28749
Immunogen	Recombinant Human Retinoblastoma-like protein 1 protein (201-500AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Retinoblastoma-like protein 1 (107 kDa retinoblastoma-associated protein) (p107) (pRb1), RBL1
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	RBL1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.