




RBM10 Antibody

Product Code	CSB-PA019400LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P98175
Immunogen	Recombinant Human RNA-binding protein 10 protein (501-700AA)
Raised In	Rabbit
Species Reactivity	Human, Mouse
Tested Applications	ELISA, WB, IHC, IF; Recommended dilution: WB:1:2000-1:10000, IHC:1:20-1:200, IF:1:50-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	RNA-binding protein 10 (G patch domain-containing protein 9) (RNA-binding motif protein 10) (RNA-binding protein S1-1) (S1-1), RBM10, DXS8237E GPATC9 GPATCH9 KIAA0122
Immunogen Species	Homo sapiens (Human)
Research Area	Epigenetics and Nuclear Signaling
Target Names	RBM10

Image



Western Blot

Positive WB detected in: Mouse heart tissue, Mouse lung tissue, Mouse brain tissue


All lanes: RBM10 antibody at 3µg/ml

Secondary


Goat polyclonal to rabbit IgG at 1/50000 dilution

Predicted band size: 104, 95, 111 kDa

Observed band size: 104 kDa



Immunohistochemistry of paraffin-embedded human lung cancer using CSB-PA019400LA01HU at dilution of 1:100



Immunofluorescent analysis of HeLa cells using CSB-PA019400LA01HU at dilution of 1:100 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)