



RBM10 Antibody, FITC conjugated

Product Code	CSB-PA019400LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P98175
Immunogen	Recombinant Human RNA-binding protein 10 protein (501-700AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	RNA-binding protein 10 (G patch domain-containing protein 9) (RNA-binding motif protein 10) (RNA-binding protein S1-1) (S1-1), RBM10, DXS8237E GPATC9 GPATCH9 KIAA0122
Immunogen Species	Homo sapiens (Human)
Research Area	Epigenetics and Nuclear Signaling
Target Names	RBM10
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.