



RBM6 Antibody, FITC conjugated

Product Code	CSB-PA019435LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P78332
Immunogen	Recombinant Human RNA-binding protein 6 protein (300-535AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Relevance	Specifically binds poly(G) RNA homopolymers in vitro.
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	RNA-binding protein 6 (Lung cancer antigen NY-LU-12) (Protein G16) (RNA-binding motif protein 6) (RNA-binding protein DEF-3), RBM6, DEF3
Species	Homo sapiens (Human)
Research Area	Epigenetics and Nuclear Signaling
Target Names	RBM6