



RCAN1 Antibody, Biotin conjugated

Product Code	CSB-PA019500LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P53805
Immunogen	Recombinant Human Calcipressin-1 protein (56-252AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Calcipressin-1 (Adapt78) (Down syndrome critical region protein 1) (Myocyte-enriched calcineurin-interacting protein 1) (MCIP1) (Regulator of calcineurin 1), RCAN1, ADAPT78 CSP1 DSC1 DSCR1
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	RCAN1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.