



RCAN1 Antibody, Biotin conjugated

Product Code	CSB-PA019500LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P53805
Immunogen	Recombinant Human Calcipressin-1 protein (56-252AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Calcipressin-1 (Adapt78) (Down syndrome critical region protein 1) (Myocyte-enriched calcineurin-interacting protein 1) (MCIP1) (Regulator of calcineurin 1), RCAN1, ADAPT78 CSP1 DSC1 DSCR1
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	RCAN1