



RCBTB2 Antibody, FITC conjugated

Product Code	CSB-PA019504LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O95199
Immunogen	Recombinant Human RCC1 and BTB domain-containing protein 2 protein (1-200AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	RCC1 and BTB domain-containing protein 2 (Chromosome condensation 1-like) (CHC1-L) (RCC1-like G exchanging factor) (Regulator of chromosome condensation and BTB domain-containing protein 2), RCBTB2, CHC1L RLG
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	RCBTB2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.