



RCBTB2 Antibody, Biotin conjugated

Product Code	CSB-PA019504LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O95199
Immunogen	Recombinant Human RCC1 and BTB domain-containing protein 2 protein (1-200AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	RCC1 and BTB domain-containing protein 2 (Chromosome condensation 1-like) (CHC1-L) (RCC1-like G exchanging factor) (Regulator of chromosome condensation and BTB domain-containing protein 2), RCBTB2, CHC1L RLG
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	RCBTB2