





REG1A Antibody

| | |
|----------------------------|---|
| Product Code | CSB-PA019546LA01HU |
| Storage | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze. |
| Uniprot No. | P05451 |
| Immunogen | Recombinant Human Lithostathine-1-alpha protein (23-166AA) |
| Raised In | Rabbit |
| Species Reactivity | Human |
| Tested Applications | ELISA, IHC; Recommended dilution: IHC:1:20-1:200 |
| Form | Liquid |
| Conjugate | Non-conjugated |
| Storage Buffer | Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4 |
| Purification Method | >95%, Protein G purified |
| Isotype | IgG |
| Clonality | Polyclonal |
| Alias | Lithostathine-1-alpha (Islet cells regeneration factor) (ICRF) (Islet of Langerhans regenerating protein) (REG) (Pancreatic stone protein) (PSP) (Pancreatic thread protein) (PTP) (Regenerating islet-derived protein 1-alpha) (REG-1-alpha) (Regenerating protein I alpha), REG1A, PSPS PSPS1 REG |
| Immunogen Species | Homo sapiens (Human) |
| Research Area | Cell Biology |
| Target Names | REG1A |

Image



Immunohistochemistry of paraffin-embedded human pancreatic tissue using CSB-PA019546LA01HU at dilution of 1:100



Immunohistochemistry of paraffin-embedded human small intestine tissue using CSB-PA019546LA01HU at dilution of 1:100