



RFPL4B Antibody, HRP conjugated

Product Code	CSB-PA019603LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q6ZWI9
Immunogen	Recombinant Human Ret finger protein-like 4B protein (1-263AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Ret finger protein-like 4B (RING finger protein 211), RFPL4B, RNF211
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	RFPL4B
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.