



RGS4 Antibody, FITC conjugated

Product Code	CSB-PA019656LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P49798
Immunogen	Recombinant Human Regulator of G-protein signaling 4 protein (1-205AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Regulator of G-protein signaling 4 (RGP4) (RGS4), RGS4
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	RGS4
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.