



ARHGEF18 Antibody, Biotin conjugated

Product Code	CSB-PA019681HD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q6ZSZ5
Immunogen	Recombinant Human Rho guanine nucleotide exchange factor 18 protein (159-551AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Rho guanine nucleotide exchange factor 18 (114 kDa Rho-specific guanine nucleotide exchange factor) (p114-Rho-GEF) (p114RhoGEF) (Septin-associated RhoGEF) (SA-RhoGEF), ARHGEF18, KIAA0521
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	ARHGEF18
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.