




RHOBTB3 Antibody

Product Code	CSB-PA019686LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O94955
Immunogen	Recombinant Human Rho-related BTB domain-containing protein 3 protein (251-500AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, WB; Recommended dilution: WB:1:1000-1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Rho-related BTB domain-containing protein 3 (EC 3.6.1.-), RHOBTB3, KIAA0878
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	RHOBTB3

Image



Western blot
All lanes: RHOBTB3 antibody at 12µg/ml + HeLa whole cell lysate
Secondary
Goat polyclonal to rabbit IgG at 1/10000 dilution
Predicted band size: 70 kDa
Observed band size: 70 kDa

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.