



RNASE9 Antibody, HRP conjugated

Product Code	CSB-PA019800LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P60153
Immunogen	Recombinant Human Inactive ribonuclease-like protein 9 protein (27-205AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Inactive ribonuclease-like protein 9, RNASE9
Immunogen Species	Homo sapiens (Human)
Research Area	Epigenetics and Nuclear Signaling
Target Names	RNASE9
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.