



RPL34 Antibody, FITC conjugated

Product Code	CSB-PA020243LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P49207
Immunogen	Recombinant Human 60S ribosomal protein L34 protein (2-117AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	60S ribosomal protein L34 (Large ribosomal subunit protein eL34), RPL34
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	RPL34
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.