



RPS26 Antibody, HRP conjugated

Product Code	CSB-PA020407LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P62854
Immunogen	Recombinant Human 40S ribosomal protein S26 protein (2-115AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	40S ribosomal protein S26 (Small ribosomal subunit protein eS26), RPS26
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	RPS26
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.