




RPS8 Antibody

Product Code	CSB-PA020478LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P62241
Immunogen	Recombinant Human 40S ribosomal protein S8 protein (2-208AA)
Raised In	Rabbit
Species Reactivity	Human, Mouse
Tested Applications	ELISA, WB; Recommended dilution: WB:1:1000-1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	40S ribosomal protein S8 (Small ribosomal subunit protein eS8), RPS8
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	RPS8

Image



Western blot

All lanes: RPS8 antibody at 6μg/ml

Lane 1: HL60 whole cell lysate

Lane 2: THP-1 whole cell lysate

Lane 3: 293T whole cell lysate

Lane 4: Mouse liver tissue

Secondary

Goat polyclonal to rabbit IgG at 1/10000 dilution

Predicted band size: 25 kDa

Observed band size: 25 kDa

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.