



RPS8 Antibody, HRP conjugated

Product Code	CSB-PA020478LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P62241
Immunogen	Recombinant Human 40S ribosomal protein S8 protein (2-208AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	40S ribosomal protein S8 (Small ribosomal subunit protein eS8), RPS8
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	RPS8