



S100A10 Antibody, HRP conjugated

Product Code	CSB-PA020623LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P60903
Immunogen	Recombinant Human Protein S100-A10 protein (2-97AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Protein S100-A10 (Calpactin I light chain) (Calpactin-1 light chain) (Cellular ligand of annexin II) (S100 calcium-binding protein A10) (p10 protein) (p11), S100A10, ANX2LG CAL1L CLP11
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	S100A10