



S100A10 Antibody, FITC conjugated

Product Code	CSB-PA020623LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P60903
Immunogen	Recombinant Human Protein S100-A10 protein (2-97AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Protein S100-A10 (Calpactin I light chain) (Calpactin-1 light chain) (Cellular ligand of annexin II) (S100 calcium-binding protein A10) (p10 protein) (p11), S100A10, ANX2LG CAL1L CLP11
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	S100A10
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.