



S100A2 Antibody, HRP conjugated

| | |
|----------------------------|---|
| Product Code | CSB-PA020630LB01HU |
| Storage | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze. |
| Uniprot No. | P29034 |
| Immunogen | Recombinant Human Protein S100-A2 protein (1-98AA) |
| Raised In | Rabbit |
| Species Reactivity | Human |
| Tested Applications | ELISA |
| Form | Liquid |
| Conjugate | HRP |
| Storage Buffer | Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4 |
| Purification Method | Antigen Affinity Purified |
| Isotype | IgG |
| Clonality | Polyclonal |
| Alias | Protein S100-A2 (CAN19) (Protein S-100L) (S100 calcium-binding protein A2), S100A2, S100L |
| Immunogen Species | Homo sapiens (Human) |
| Research Area | Signal Transduction |
| Target Names | S100A2 |