



S100A2 Antibody, Biotin conjugated

Product Code	CSB-PA020630LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P29034
Immunogen	Recombinant Human Protein S100-A2 protein (1-98AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Protein S100-A2 (CAN19) (Protein S-100L) (S100 calcium-binding protein A2), S100A2, S100L
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	S100A2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.