



SAG Antibody, HRP conjugated

Product Code	CSB-PA020669YB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P10523
Immunogen	Recombinant Human S-arrestin protein (1-405AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	S-arrestin (48 kDa protein) (Retinal S-antigen) (S-AG) (Rod photoreceptor arrestin), SAG
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	SAG
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.