



SCAMP2 Antibody, Biotin conjugated

Product Code	CSB-PA020751LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O15127
Immunogen	Recombinant Human Secretory carrier-associated membrane protein 2 protein (1-153AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Secretory carrier-associated membrane protein 2 (Secretory carrier membrane protein 2), SCAMP2
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	SCAMP2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.