



# SCGB2A1 Antibody, HRP conjugated

Product Code	CSB-PA020816LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O75556
Immunogen	Recombinant Human Mammaglobin-B protein (19-95AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Mammaglobin-B (Lacryglobin) (Lipophilin-C) (Mammaglobin-2) (Secretoglobin family 2A member 1), SCGB2A1, LIPHC MGB2 UGB3
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	SCGB2A1